

**Respiratory Pathogen Testing?  
So many SWABS!  
START HERE**

BBL culture swab  
(double red cap)  
PS# 0500356



← Group A Strep

**WHAT AM I TESTING FOR?**

Respiratory Pathogen Panel (RPP) and/or Covid Molecular

BD Universal Viral Transport  
PS# 5017049



Respiratory Syncytial Virus (RSV)  
or  
Influenza (FLU): SARS Covid Rapid Antigen

**ARE YOU PERFORMING POINT OF CARE OR SENDING IT TO A LAB?**

POC

Sending to a lab

**FLU/SARS OR RSV**

**WHERE ARE YOU LOCATED?**

FLU/SARS

RSV


WH, UH or clinic with lab services

MOI or clinic without lab service


Use the foam tipped nasal swab  
**Included in kit**





Use the Hydraflock Sterile flocked nasopharyngeal swab  
**Included in kit**



**Option 1 for FLU/SARS Combo:**  
Use one foam tipped swab from kit  
**Option 2 for FLU/SARS combo:** use one copan FLOQ N.P. swab  
**Options for RSV:** use one swab from kit or use one copan FLOQ N.P. swab PS# 5027950  
**Place in dry tube and send to lab ASAP**



Use a FLOQ nasopharyngeal swab.  
Place in 1 ml saline in sterile tube  
Swab PS# 5019619  
Sterile Tube PS# 5007415  
Saline PS# 5003758

Always refer to the Pathology Test Catalog at [muhealth.testcatalog.org](http://muhealth.testcatalog.org) found on Citrix receiver or in Power Chart for the most update-to-date testing information.

*For convenience, e-procure numbers have been added to this poster for ordering purposes.*